











Vessel Format	Type of tips	epT.I.P.S.®			epT.I.P.S.® 384	
	Electronic Pipette	Eppendorf Xplorer® plus Move It®			Eppendorf Xplorer® plus Move It®	
	Mechanical Pipette	Eppendorf Research® plus Move It®			Eppendorf Research® plus Move It®	
	No. of channels	4	6	8	8	12
	Volume (µL)	300 / 1,200	300 / 1,200	300 / 1,200	20 / 100	20 / 100
	Tip distance (mm)	9-33	9-20	9-14	4.5-14	4.5-9
384 Wells 	4.5				■	■
96 Wells / PCR Tube Strips 	9	■	■	■	■	■
48 Wells 	13	■	■	■	■	
24 Wells 	19	■	■			
12 Wells 	26	■				
1.5 / 2.0 / 5.0 / 15 mL Tubes 	21-33	■				
1.5 / 2.0 / 5.0 / 15 mL Tubes 	15-20	■	■			
0.5 / 1.5 / 2.0 mL Tubes 	9-14	■	■	■	■	
Agarose gel 	~ 4.5-9	■ *	■ *	■ *	■	■

* Limited suitability due to volume and size of tips